



## CytoBrite® *PLUS* Slide Incubation System

Only FISH Hybridizer Fully Compliant with New CAP Regulation CYG.33950



Designed to meet the requirements of new CAP regulation CYG.33950, the CytoBrite® *PLUS* Slide Incubation System is equipped with a separate temperature sensor for each of its 12 slide positions. Simply connect the provided NIST-certified digital thermometer to view each slide position temperature.



Cytogenetics Checklist | 08.17.2016

**\*\*NEW\*\* CAP REGULATION**

**CYG.33950**

**ISH Slide Processing System Temperature Checks**

**Phase II**

Individual slide slots (or a representative sample thereof) of *in situ* hybridization (ISH) temperature controlled slide processing systems are checked for temperature accuracy before being placed in service and at least annually thereafter.

**Evidence of compliance:**

- √ Written procedure for verification of temperature accuracy **AND**
- √ Records of equipment verification

