

## Model 2000 Hybridization Incubator

*Small and Economical for Lower Volume Usage*

The Model 2000 Hybridization Incubator is designed for laboratories that process a limited number of membranes, have limited space or budget or require multiple units operating at different temperatures.



Ideal for crowded laboratories, the compact Model 2000 Hybridization Incubator occupies less than 1.5 sq. ft. of bench space and can be stacked for additional space savings. The Model 2000 uses a high energy heating element with a fan that brings the chamber to the desired temperature in minutes. The precision heat controller effectively regulates the temperature to  $\pm 0.10^{\circ}\text{C}$  and can be toggled to display either the set or actual temperature. The rotator when activated moves at a fixed 6 RPM and is controlled with a three position switch (on, off, jog) that can be used to momentarily advance the rotator to add or remove tubes. The standard rotator provided with the unit holds up to four, long (300mm) or eight, short (150mm) glass tubes. One long (300mm) glass tube with cap is provided with the system. Additional tubes and racks are available separately. Priced at less than half of larger models, the Model 2000 is ideal for the budget minded laboratory.

### Model 2000 Hybridization Incubator Specifications

<b>Tube Capacity</b>	4–35x300mm 8–35x150mm
<b>Temperature Range/ Precision</b>	Ambient +10 to 99°C ( $\pm 0.10^{\circ}\text{C}$ )
<b>Rotator Speed</b>	6 rpm
<b>Dimensions</b>	Exterior: 13.0x16.0x14 inches (33.0x40.5x35.5cm) Chamber: 8.0x14.0x7.5 inches (20.0x35.5x35.5cm)
<b>Weight (without Rotator)</b>	32 lbs. (14.5kg)0

### Ordering Information

Model 2000 Hybridization Incubator and Accessories		
Part No.	Description	Unit
1040-50-1	Model 2000 Hybridization Incubator, 115V. Includes standard 4/8 tube rotator and one 35x300mm tube with cap.	Each
1040-50-2	Model 2000 Hybridization Incubator, 220V. Includes standard 4/8 tube rotator and one 35x300mm tube with cap.	Each
1040-51-1	Standard 4/8 tube rotator for Model 2000 Hybridization Incubator.	Each
1040-20-1	Replacement drip pan for Model 2000 Hybridization Incubator.	Each
Hybridization Tubes and Other Incubator Supplies		
Part No.	Description	Unit
1040-01-0	Hybridization tube, 35x300mm with screw-cap.	Each
1040-02-0	Hybridization tube, 35x150mm with screw-cap.	Each
1040-05-0	Hybridization tube, 35x300mm without screw-cap.	Each
1040-06-0	Hybridization tube, 35x150mm without screw-cap.	Each
1040-65-0	Tube gripper, textured rubber pad	10/Pack
1051-52-0	Handheld electronic digital thermometer with Type T thermocoupler—includes NIST traceable certificate of calibration and interface cable	Each

### Product Sales and Support Contact Information

#### Mailing Address:

SciGene  
306 Potrero Avenue  
Sunnyvale, CA 94085  
USA

Business Hours: 8:00 AM–5:00 PM Pacific Time

Phone: 408-733-7337/800-342-2119

Fax: 408-733-7336

#### e-mail:

Customer Service: [custserv@scigene.com](mailto:custserv@scigene.com)

Technical Service: [techserv@scigene.com](mailto:techserv@scigene.com)

Sales Information: [sales@scigene.com](mailto:sales@scigene.com)

Website: [www.scigene.com](http://www.scigene.com)